# Metabolite annotation by Liquid-Chromatography-Mass Spectrometry (LC-MS) Practice

#### 1. Overview

This practice emulates the core functionality of the metabolite annotation tools, which integrates liquid chomatography and mass spectrometry data with metabolite databases to identify potential compounds and adducts based on the data acquired by analytical instrumentation. Concretely it will be based on the Retention Time (RT) acquired by Liquid Chromatography (LC) and the experimental mass-to-charge ratio (m/z) obtained by Mass Spectrometry (MS) means. The template code will be a J2EE project using:

- \*\*Drools\*\*: For rule-based reasoning to match m/z values to metabolite candidates and subsequently use the elution order based on the RT to confirm or reject annotations.
- \*\*J2EE\*\*: For application structure.
- \*\*Maven\*\*: For dependency management and build automation.
- \*\*JUnit + Mockito\*\*: For unit testing.

### 2. Architecture

mass-spec-ceu-practice/

```
| ├─ main/

| | ├─ java/

| | | ├─ resources/

| | ├─ rules/

| | ├─ metabolite-matching.drl

| └─ test/

| └─ java


| └─ lipid/AdductDetectionTest.java

| └─ lipid/ElutionOrderTest.java

| ├─ pom.xml
```

## 3. Sample Model

The proposed UML is simplified to focus



Where the class Annotation is the data type which will be used as a fact. It contains the next attributes

- intensity: double |
- rtMin: double |

```
| - adduct: String
| - groupedSignals: Set<Peak>|
- score: int
+ Annotation(...)
| + getLipid(): Lipid
| + getMz(): double
| + getRtMin(): double
| + getAdduct(): String
| + setAdduct(String): void |
| + getIntensity(): double |
| + getGroupedSignals(): Set<Peak> |
| + getScore(): int
| + setScore(int): void |
| + addScore(int): void |
| + equals(Object): boolean |
+ hashCode(): int
| + toString(): String
```

Check the project template and go to the Junit test classes. There are 2 main goals:

1- Adduct detection. Based on the masses difference between the grouped mzs obtained in a corresponding RT. The mass differences will be used to detect the adduct based on the mass difference within a certain tolerance specified in mDa or in ppm. The ppm

```
@Test
public void shouldDetectAdductBasedOnMzDifference() {

// Given two peaks with ~21.98 Da difference (e.g., [M+H]+ and [M+Na]+)
Peak mH = new Peak(700.500, 1000000.0); // [M+H]+
Peak mNa = new Peak(722.482, 80000.0); // [M+Na]+
```

```
Lipid lipid = new Lipid(1, "PC 34:1", "C42H82NO8P", "PC", 34, 1);
  double annotationMZ = 700.49999d;
  double annotationIntensity = 80000.0;
  double annotationRT = 6.5d;
  Annotation annotation = new Annotation(lipid, annotationMZ, annotationIntensity,
annotationRT, Set.of(mH, mNa));
  // Then we should call the algorithmic/knowledge system rules fired to detect the
adduct and Set it!
  assertNotNull("[M+H]+ should be detected", annotation.getAdduct());
  assertEquals( "Adduct inferred from lowest mz in group", "[M+H]+",
annotation.getAdduct());
@Test
public void shouldDetectLossOfWaterAdduct() {
                                             // [M+H]+
  Peak mh = new Peak(700.500, 90000.0);
  Peak mhH2O = new Peak(682.4894, 70000.0); // [M+H-H_2O]+, ~18.0106 Da
less
  Lipid lipid = new Lipid(1, "PE 36:2", "C41H78NO8P", "PE", 36, 2);
  Annotation annotation = new Annotation(lipid, mh.getMz(), mh.getIntensity(),
7.5d, Set.of(mh, mhH2O));
  assertNotNull("[M+H]+ should be detected", annotation.getAdduct());
  assertEquals( "Adduct inferred from lowest mz in group", "[M+H]+",
annotation.getAdduct());
@Test
public void shouldDetectDoublyChargedAdduct() {
  // Assume real M = (700.500 - 1.0073) = 699.4927
  // So [M+2H]2+ = (M + 2.0146) / 2 = 350.7536
  Peak singlyCharged = new Peak(700.500, 100000.0); // [M+H]+
  Peak doublyCharged = new Peak(350.754, 85000.0); // [M+2H]2+
  Lipid lipid = new Lipid(3, "TG 54:3", "C57H104O6", "TG", 54, 3);
  Annotation annotation = new Annotation(lipid, singlyCharged.getMz(),
singlyCharged.getIntensity(), 10d, Set.of(singlyCharged, doublyCharged));
  assertNotNull("[M+H]+ should be detected", annotation.getAdduct());
```

```
assertEquals( "Adduct inferred from lowest mz in group","[M+H]+",
annotation.getAdduct());
}
```

#### Elution order of compounds:

```
@Test
public void score1BasedOnRT() {
  // Assume lipids already annotated
  LOG.info("Creating RuleUnit");
  LipidScoreUnit lipidScoreUnit = new LipidScoreUnit();
  RuleUnitInstance<LipidScoreUnit> instance =
RuleUnitProvider.get().createRuleUnitInstance(lipidScoreUnit);
  // TODO CHECK THE Monoisotopic MASSES OF THE COMPOUNDS IN
https://chemcalc.org/
  Lipid lipid1 = new Lipid(1, "TG 54:3", "C57H104O6", "TG", 54, 3); // MZ of [M+H]+
= 885.79057
  Lipid lipid2 = new Lipid(2, "TG 52:3", "C55H100O6", "TG", 52, 3); // MZ of [M+H]+
= 857.75927
  Lipid lipid3 = new Lipid(3, "TG 56:3", "C59H108O6", "TG", 56, 3); // MZ of [M+H]+
= 913.82187
  Annotation annotation1 = new Annotation(lipid1, 885.79056, 10E6, 10d);
  Annotation annotation2 = new Annotation(lipid2, 857.7593, 10E7, 9d);
  Annotation annotation3 = new Annotation(lipid3, 913.822, 10E5, 11d);
  LOG.info("Insert data");
  try {
    lipidScoreUnit.getAnnotations().add(annotation1);
    lipidScoreUnit.getAnnotations().add(annotation2);
    lipidScoreUnit.getAnnotations().add(annotation3);
    LOG.info("Run query. Rules are also fired");
    instance.fire();
    // Here the logic that we expect. In this case we expect the full 3 annotations to
have a positive score of 1
    assertEquals(1.0, annotation1.getNormalizedScore(), 0.01);
    assertEquals(1.0, annotation2.getNormalizedScore(), 0.01);
    assertEquals(1.0, annotation3.getNormalizedScore(), 0.01);
```

```
}
finally {
  instance.close();
}
```

// TODO NEW RULES WILL BE ADDED HERE FOR THE RT-ELUTION ORDER. The rules are defined in the section 4.

## 4. Drools Rules (metabolite-matching.drl)

This is a dummy example of a rule. This rule is the basis of the RT rules, but the rules will be defined after so you can translate them into your project.

```
// TODO Include here rules and queries to fulfill the practice requirements
// This is one example of rules that only prints the factorial combination of all pairs
eliminating the A-A combination.
rule "Score 1 for lipid pair with increasing RT and carbon count"
when
  $a1 : /annotations [$rt1 : rtMin, $carbonCount1 : lipid.getCarbonCount(),
$doubleBondCount : lipid.getDoubleBondsCount()]
  $a2 : /annotations [this!= $a1, lipid.getDoubleBondsCount() ==
$doubleBondCount, rtMin > $rt1, lipid.getCarbonCount() > $carbonCount1]
  // in this case, the only change is the addition of the score, but the fact does not
change so we do not break the principle of refractoriness
  $a1.addScore(1);
  $a2.addScore(1);
  //!! TODO ONLY FOR DEBUGGING
  System.out.println("Scored 1 for annotations: " + $a1 + " and " + $a2);
end
```

Here the rules to fulfill the Tests stated in section 8 should be fulfilled:

- 1- Adduct detection based on the mass differences of the grouped peaks
- 2- RT Elution order based on positive evidence. The score should be positive according to the methods "addScore(int delta)" and "getNormalizedScore()" of the class Annotation.

Each rule adds a score of 1.

- 1. RT(L1) > RT (L2) if L1 Lipid Type == L2 LypidType; L1 Double Bonds == L2 Double bonds and L1 carbons > L2 carbons.
- 2. RT(L1) > RT (L2) if L1 Lipid Type == L2 LypidType; L1 Double Bonds < L2 Double bonds and L1 carbons == L2 carbons.
- RT(L1) > RT (L2) if L1 Lipid Type != L2 LypidType; L1 Double Bonds == L2 Double bonds and L1 carbons == L2 carbons according to this pattern: -> PG < PE < PI < PA < PS << PC.</li>
- 2- RT Elution order based on negative evidence.

Each rule add a score of -1.

- 4. RT(L1) < RT (L2) if L1 Lipid Type == L2 LypidType; L1 Double Bonds == L2 Double bonds and L1 carbons > L2 carbons.
- 5. RT(L1) < RT (L2) if L1 Lipid Type == L2 LypidType; L1 Double Bonds < L2 Double bonds and L1 carbons == L2 carbons.
- 6. RT(L1) < RT (L2) if L1 Lipid Type != L2 LypidType; L1 Double Bonds == L2 Double bonds and L1 carbons == L2 carbons according to this pattern: -> PG < PE < PI < PA < PS << PC.

## 5. Rule Engine Service

The student can choose among the classical drools KieSessions or the ruleUnitData. The example of the practice currently is performed with RuleUnitData to keep the consistency with the previous examples.

This is an example of the KieSession instead of the RuleUnitData (Check the templatePractice or all the examples provided in this course).

package droolsService;

import Lipid;

import org.kie.api.KieServices;

import org.kie.api.runtime.KieContainer;

import org.kie.api.runtime.KieSession;

import java.util.List;

```
public class RuleEngineService {
    private final KieContainer kieContainer;

public RuleEngineService() {
    this.kieContainer = KieServices.Factory.get().getKieClasspathContainer();
}

public void applyRules(List<Lipid> lipids) {
    KieSession kieSession = kieContainer.newKieSession();
    for (Lipid l: lipids) {
        kieSession.insert(l);
    }
    kieSession.fireAllRules();
    kieSession.dispose();
}
```

Here is the example with RuleUnitData.

RuleUnitData example:

```
LipidScoreUnit lipidScoreUnit = new LipidScoreUnit();

RuleUnitInstance < LipidScoreUnit > instance = RuleUnitProvider.get().createRuleUnitInstance(lipidScoreUnit);

Lipid lipid1 = new Lipid(1, "TG 54:3", "C57H104O6", "TG", 54, 3); // MZ of [M+H]+ = 885.79057

Lipid lipid2 = new Lipid(2, "TG 52:3", "C55H100O6", "TG", 52, 3); // MZ of [M+H]+ = 857.75927

Lipid lipid3 = new Lipid(3, "TG 56:3", "C59H108O6", "TG", 56, 3); // MZ of [M+H]+ =
```

```
913.82187
Annotation annotation1 = new Annotation(lipid1, 885.79056, 10E6, 10d);
Annotation annotation2 = new Annotation(lipid2, 857.7593, 10E7, 9d);
Annotation annotation3 = new Annotation(lipid3, 913.822, 10E5, 11d);

LOG.info("Insert data");

try {
    lipidScoreUnit.getAnnotations().add(annotation1);
    lipidScoreUnit.getAnnotations().add(annotation2);
    lipidScoreUnit.getAnnotations().add(annotation3);

LOG.info("Run query. Rules are also fired");
    instance.fire();

}
finally {
    instance.close();
}
```

## 7. Maven Dependencies (pom.xml)

```
<?xml version="1.0" encoding="UTF-8"?>
<project xmlns="http://maven.apache.org/POM/4.0.0"
xmlns:xsi="http://www.w3.org/2001/XMLSchema-instance"
xsi:schemaLocation="http://maven.apache.org/POM/4.0.0
http://maven.apache.org/xsd/maven-4.0.0.xsd">
<modelVersion>4.0.0</modelVersion>

<groupId>dss</groupId>
<artifactId>template</artifactId>
<version>1.0-SNAPSHOT</version>
<packaging>kjar</packaging>

<name>template</name>

<
```

```
<drools-version>10.0.0</drools-version>
 <junit-version>4.13.2</junit-version>
 logback-version>1.5.17</logback-version>
 <slf4j-version>2.0.17</slf4j-version>
</properties>
<dependencyManagement>
<dependencies>
  <dependency>
   <groupId>org.drools
   <artifactId>drools-bom</artifactId>
   <type>pom</type>
   <version>${drools-version}
  </dependency>
 </dependencies>
</dependencyManagement>
<dependencies>
<dependency>
  <groupId>org.drools</groupId>
  <artifactId>drools-ruleunits-engine</artifactId>
  <version>${drools-version}
 </dependency>
 <dependency>
  <groupId>junit</groupId>
  <artifactId>junit</artifactId>
  <version>${junit-version}</version>
  <scope>test</scope>
 </dependency>
 <dependency>
  <groupId>org.slf4j</groupId>
  <artifactId>slf4j-api</artifactId>
  <version>${slf4j-version}
 </dependency>
 <dependency>
  <groupId>ch.qos.logback
  <artifactId>logback-classic</artifactId>
  <version>${logback-version}</version>
 </dependency>
 <dependency>
```

```
<groupId>org.drools</groupId>
  <artifactId>drools-wiring-dynamic</artifactId>
  <version>${drools-version}</version>
 </dependency>
</dependencies>
<bul>duild>
 <plugins>
  <plugin>
   <artifactId>maven-compiler-plugin</artifactId>
   <version>${maven-compiler-version}</version>
   <configuration>
    <release>${maven.compiler.release}</release>
   </configuration>
  </plugin>
  <plugin>
   <groupId>org.kie</groupId>
   <artifactId>kie-maven-plugin</artifactId>
   <version>${drools-version}
   <extensions>true</extensions>
  </plugin>
 </plugins>
</build>
</project>
```

#### 8. Unit Tests

Find them in the template project. These unit tests contain the requirements of the practice. There are two type of knowledge here:

REMINDER FROM Section 4.

- 1- Adduct detection based on the mass differences of the grouped peaks
- 2- RT Elution order based on positive evidence. The score should be positive according to the methods "addScore(int delta)" and "getNormalizedScore()" of the class Annotation.

Each rule adds a score of 1.

7. RT(L1) > RT (L2) if L1 Lipid Type == L2 LypidType; L1 Double Bonds == L2 Double bonds and L1 carbons > L2 carbons.

- 8. RT(L1) > RT (L2) if L1 Lipid Type == L2 LypidType; L1 Double Bonds < L2 Double bonds and L1 carbons == L2 carbons.
- RT(L1) > RT (L2) if L1 Lipid Type != L2 LypidType; L1 Double Bonds == L2 Double bonds and L1 carbons == L2 carbons according to this pattern: -> PG < PE < PI < PA < PS << PC.</li>
- 2- RT Elution order based on negative evidence.

Each rule add a score of -1.

- 10. RT(L1) < RT (L2) if L1 Lipid Type == L2 LypidType; L1 Double Bonds == L2 Double bonds and L1 carbons > L2 carbons.
- 11. RT(L1) < RT (L2) if L1 Lipid Type == L2 LypidType; L1 Double Bonds < L2 Double bonds and L1 carbons == L2 carbons.
- 12. RT(L1) < RT (L2) if L1 Lipid Type != L2 LypidType; L1 Double Bonds == L2 Double bonds and L1 carbons == L2 carbons according to this pattern: -> PG < PE < PI < PA < PS << PC.

Those are the Junit Tests to pass. BlackBox Unit Tests will be run to check the proper functioning of the application.

```
package lipid;
import org.drools.ruleunits.api.RuleUnitInstance;
import org.drools.ruleunits.api.RuleUnitProvider;
import org.junit.Before;
import org.junit.Test;
import org.slf4j.Logger;
import org.slf4j.LoggerFactory;
import java.util.List;
import java.util.Set;
import static org.junit.Assert.assertEquals;
import static org.junit.Assert.assertNotNull;
import static org.junit.Assert.assertTrue;
public class ElutionOrderTest {

static final Logger LOG = LoggerFactory.getLogger(ElutionOrderTest.class);
// !!TODO For the adduct detection both regular algorithms or drools can be used
```

```
as far the tests are passed.
  @Before
  public void setup() {
    //!! TODO Empty by now, you can create common objects for all tests.
  * Test to check the elution order of the lipids. The elution order is based on the
number of carbons if the lipid type and the number of
   st double bonds is the same. The larger the number of carbons, the longer the RT.
  @Test
  public void score1BasedOnRTCarbonNumbers() {
    LOG.info("Creating RuleUnit");
    LipidScoreUnit lipidScoreUnit = new LipidScoreUnit();
    RuleUnitInstance<LipidScoreUnit> instance =
RuleUnitProvider.get().createRuleUnitInstance(lipidScoreUnit);
    // TODO CHECK THE Monoisotopic MASSES OF THE COMPOUNDS IN
https://chemcalc.org/
    Lipid lipid1 = new Lipid(1, "TG 54:3", "C57H104O6", "TG", 54, 3); // MZ of [M+H]
+ = 885.79057
    Lipid lipid2 = new Lipid(2, "TG 52:3", "C55H100O6", "TG", 52, 3); // MZ of [M+H]
+ = 857.75927
    Lipid lipid3 = new Lipid(3, "TG 56:3", "C59H108O6", "TG", 56, 3); // MZ of [M+H]
+ = 913.82187
    Annotation annotation1 = new Annotation(lipid1, 885.79056, 10E6, 10d,
loniationMode.POSITIVE);
    Annotation annotation2 = new Annotation(lipid2, 857.7593, 10E7, 9d,
loniationMode.POSITIVE);
    Annotation annotation3 = new Annotation(lipid3, 913.822, 10E5, 11d,
loniationMode.POSITIVE);
    LOG.info("Insert data");
    try {
       lipidScoreUnit.getAnnotations().add(annotation1);
       lipidScoreUnit.getAnnotations().add(annotation2);
       lipidScoreUnit.getAnnotations().add(annotation3);
       LOG.info("Run query. Rules are also fired");
       instance.fire():
```

```
// Here the logic that we expect. In this case we expect the full 3 annotations
to have a positive score of 1
       assertEquals(1.0, annotation1.getNormalizedScore(), 0.01);
       assertEquals(1.0, annotation2.getNormalizedScore(), 0.01);
       assertEquals(1.0, annotation3.getNormalizedScore(), 0.01);
    finally {
       instance.close();
   *!!TODO Test to check the elution order of the lipids. The elution order is based
on the number of double bonds if the lipid type and the number of
  @Test
  public void score1BasedOnRTDoubleBonds() {
    // Assume lipids already annotated
    LOG.info("Creating RuleUnit");
    LipidScoreUnit lipidScoreUnit = new LipidScoreUnit();
    RuleUnitInstance<LipidScoreUnit> instance =
RuleUnitProvider.get().createRuleUnitInstance(lipidScoreUnit);
    // TODO CHECK THE Monoisotopic MASSES OF THE COMPOUNDS IN
https://chemcalc.org/
    Lipid lipid1 = new Lipid(1, "TG 54:3", "C57H104O6", "TG", 54, 3); // MZ of [M+H]
+ = 885.79057
    Lipid lipid2 = new Lipid(2, "TG 54:4", "C57H102O6", "TG", 54, 4); // MZ of [M+H]
+ = 883.77492
    Lipid lipid3 = new Lipid(3, "TG 54:2", "C57H106O6", "TG", 54, 2); // MZ of [M+H]
+ = 887.80622
    Annotation annotation1 = new Annotation(lipid1, 885.79056, 10E6, 10d,
IoniationMode.POSITIVE);
     Annotation annotation2 = new Annotation(lipid2, 883.77492, 10E7, 9d,
loniationMode.POSITIVE);
    Annotation annotation3 = new Annotation(lipid3, 887.80622, 10E5, 11d,
IoniationMode.POSITIVE);
    LOG.info("Insert data");
```

```
lipidScoreUnit.getAnnotations().add(annotation1);
       lipidScoreUnit.getAnnotations().add(annotation2);
       lipidScoreUnit.getAnnotations().add(annotation3);
       LOG.info("Run query. Rules are also fired");
       instance.fire();
       // Here the logic that we expect. In this case we expect the full 3 annotations
to have a positive score of 1
       assertEquals(1.0, annotation1.getNormalizedScore(), 0.01);
       assertEquals(1.0, annotation2.getNormalizedScore(), 0.01);
       assertEquals(1.0, annotation3.getNormalizedScore(), 0.01);
    finally {
       instance.close();
  * !!TODO Test to check the elution order of the lipids. The elution order is based
on the number of double bonds if the lipid type and the number of
   st carbons is the same. The higher the number of double bonds, the shorter the RT.
  * The RT order of lipids with the same number of carbons and double bonds is the
same
   * -> PG < PE < PI < PA < PS << PC.
  @Test
  public void score1BasedOnRTDoubleBonds() {
    // Assume lipids already annotated
    LOG.info("Creating RuleUnit");
    LipidScoreUnit lipidScoreUnit = new LipidScoreUnit();
    RuleUnitInstance<LipidScoreUnit> instance =
RuleUnitProvider.get().createRuleUnitInstance(lipidScoreUnit);
    // TODO CHECK THE Monoisotopic MASSES OF THE COMPOUNDS IN
https://chemcalc.org/
    Lipid lipid1 = new Lipid(1, "PI 34:0", "C43H83O13P", "PI", 54, 0); // MZ of [M+H]
+ = 839.56441
     Lipid lipid2 = new Lipid(2, "PG 34:0", "C40H79O10P", "PG", 54, 0); // MZ of
[M+H]+ = 751.54836
     Lipid lipid3 = new Lipid(3, "PC 34:0", "C42H84NO8P", "PC", 54, 0); // MZ of
[M+H]+ = 762.60073
     Annotation annotation 1 = \text{new Annotation}(\text{lipid1}, 839.5644179056, 10E6, 10d,
```

```
Annotation annotation2 = new Annotation(lipid2, 751.54836, 10E7, 9d,
IoniationMode.POSITIVE);
     Annotation annotation3 = new Annotation(lipid3, 913.822, 10E5, 11d,
IoniationMode.POSITIVE):
    LOG.info("Insert data");
    try {
       lipidScoreUnit.getAnnotations().add(annotation1);
       lipidScoreUnit.getAnnotations().add(annotation2);
       lipidScoreUnit.getAnnotations().add(annotation3);
       LOG.info("Run guery. Rules are also fired");
       instance.fire();
       // Here the logic that we expect. In this case we expect the full 3 annotations
to have a positive score of 1
       assertEquals(1.0, annotation1.getNormalizedScore(), 0.01);
       assertEquals(1.0, annotation2.getNormalizedScore(), 0.01);
       assertEquals(1.0, annotation3.getNormalizedScore(), 0.01);
    finally {
       instance.close();
  //!! TODO. PART II OF THE PRACTICE
  //!! TODO FROM HERE. This part of the practice is not mandatory to pass the
course. It provides negative evidence when the RT elution order does not fulfill the
patterns.
  *!!TODO Test to check the elution order of the lipids. The elution order is based
on the number of double bonds if the lipid type and the number of
   st carbons is the same. The higher the number of double bonds, the shorter the RT.
same
   * -> PG < PE < PI < PA < PS << PC.
  @Test
  public void score1BasedOnRTDoubleBonds() {
    // Assume lipids already annotated
    LOG.info("Creating RuleUnit");
    LipidScoreUnit lipidScoreUnit = new LipidScoreUnit();
    RuleUnitInstance<LipidScoreUnit> instance =
```

```
RuleUnitProvider.get().createRuleUnitInstance(lipidScoreUnit);
    // TODO CHECK THE Monoisotopic MASSES OF THE COMPOUNDS IN
https://chemcalc.org/
     Lipid lipid1 = new Lipid(1, "PI 34:0", "C43H83O13P", "PI", 54, 0); // MZ of [M+H]
+ = 839.56441
    Lipid lipid2 = new Lipid(2, "PG 34:0", "C40H79O10P", "PG", 54, 0); // MZ of
[M+H]+ = 751.54836
    Lipid lipid3 = new Lipid(3, "PC 34:0", "C42H84N08P", "PC", 54, 0); // MZ of
[M+H]+ = 762.60073
     Annotation annotation1 = new Annotation(lipid1, 839.5644179056, 10E6, 10d,
IoniationMode.POSITIVE):
    Annotation annotation2 = new Annotation(lipid2, 751.54836, 10E7, 9d,
loniationMode.POSITIVE);
     Annotation annotation3 = new Annotation(lipid3, 913.822, 10E5, 8d,
loniationMode.POSITIVE);
    LOG.info("Insert data");
    try {
       lipidScoreUnit.getAnnotations().add(annotation1);
       lipidScoreUnit.getAnnotations().add(annotation2);
       lipidScoreUnit.getAnnotations().add(annotation3);
       LOG.info("Run query. Rules are also fired");
       instance.fire();
       // Here the logic that we expect. In this case we expect the full 3 annotations
to have a positive score of 1
       assertEquals(0d, annotation1.getNormalizedScore(), 0.01); // !! !! TODO the
scores should be between -1 and 1. It is done, but check it out for yourself
       assertEquals(0d, annotation2.getNormalizedScore(), 0.01); // !! TODO the
scores should be between -1 and 1. It is done, but check it out for yourself
       assertEquals(-1.0, annotation3.getNormalizedScore(), 0.01); // !! TODO the
scores should be between -1 and 1. It is done, but check it out for yourself
    finally {
       instance.close();
  *!!TODO Test to check the elution order of the lipids. The elution order is based
on the number of double bonds if the lipid type and the number of
```

```
st carbons is the same. The higher the number of double bonds, the shorter the RT.
  @Test
  public void score1BasedOnRTDoubleBonds() {
    // Assume lipids already annotated
    LOG.info("Creating RuleUnit");
    LipidScoreUnit lipidScoreUnit = new LipidScoreUnit();
    RuleUnitInstance<LipidScoreUnit> instance =
RuleUnitProvider.get().createRuleUnitInstance(lipidScoreUnit);
    // TODO CHECK THE Monoisotopic MASSES OF THE COMPOUNDS IN
https://chemcalc.org/
    Lipid lipid1 = new Lipid(1, "TG 54:3", "C57H104O6", "TG", 54, 3); // MZ of [M+H]
+ = 885.79057
    Lipid lipid2 = new Lipid(2, "TG 54:4", "C57H102O6", "TG", 54, 4); // MZ of [M+H]
    Lipid lipid3 = new Lipid(3, "TG 54:2", "C57H106O6", "TG", 54, 2); // MZ of [M+H]
+ = 887.80622
     Annotation annotation 1 = \text{new Annotation}(\text{lipid1}, 885.79056, 10E6, 10d,
loniationMode.POSITIVE);
    Annotation annotation2 = new Annotation(lipid2, 883.77492, 10E7, 9d,
loniationMode.POSITIVE);
    Annotation annotation3 = new Annotation(lipid3, 887.80622, 10E5, 8d,
loniationMode.POSITIVE);
    LOG.info("Insert data");
    try {
       lipidScoreUnit.getAnnotations().add(annotation1);
       lipidScoreUnit.getAnnotations().add(annotation2);
       lipidScoreUnit.getAnnotations().add(annotation3);
       LOG.info("Run query. Rules are also fired");
       instance.fire();
       // Here the logic that we expect. In this case we expect the full 3 annotations
to have a positive score of 1
       assertEquals(0d, annotation1.getNormalizedScore(), 0.01); // !! TODO the
scores should be between -1 and 1. It is done, but check it out for yourself
       assertEquals(0d, annotation2.getNormalizedScore(), 0.01); // !! TODO the
scores should be between -1 and 1. It is done, but check it out for yourself
       assertEquals(-1.0, annotation3.getNormalizedScore(), 0.01); // !! !! TODO the
scores should be between -1 and 1. It is done, but check it out for yourself
```

```
finally {
       instance.close();
   * Test to check the elution order of the lipids. The elution order is based on the
number of carbons if the lipid type and the number of
   st double bonds is the same. The larger the number of carbons, the longer the RT.
  @Test
  public void score1BasedOnRTCarbonNumbers() {
    // Assume lipids already annotated
    LOG.info("Creating RuleUnit");
    LipidScoreUnit lipidScoreUnit = new LipidScoreUnit();
     RuleUnitInstance < LipidScoreUnit > instance =
RuleUnitProvider.get().createRuleUnitInstance(lipidScoreUnit);
     // TODO CHECK THE Monoisotopic MASSES OF THE COMPOUNDS IN
https://chemcalc.org/
    Lipid lipid1 = new Lipid(1, "TG 54:3", "C57H104O6", "TG", 54, 3); // MZ of [M+H]
+ = 885.79057
    Lipid lipid2 = new Lipid(2, "TG 52:3", "C55H100O6", "TG", 52, 3); // MZ of [M+H]
+ = 857.75927
    Lipid lipid3 = new Lipid(3, "TG 56:3", "C59H108O6", "TG", 56, 3); // MZ of [M+H]
+ = 913.82187
    Annotation annotation1 = new Annotation(lipid1, 885.79056, 10E6, 10d,
IoniationMode.POSITIVE);
     Annotation annotation2 = new Annotation(lipid2, 857.7593, 10E7, 9d,
IoniationMode.POSITIVE):
    Annotation annotation3 = new Annotation(lipid3, 913.822, 10E5, 8d,
loniationMode.POSITIVE);
    LOG.info("Insert data");
    try {
       lipidScoreUnit.getAnnotations().add(annotation1);
       lipidScoreUnit.getAnnotations().add(annotation2);
       lipidScoreUnit.getAnnotations().add(annotation3);
       LOG.info("Run query. Rules are also fired");
       instance.fire();
       // Here the logic that we expect. In this case we expect the full 3 annotations
to have a positive score of 1
```

```
assertEquals(0d, annotation1.getNormalizedScore(), 0.01); // !! TODO the
scores should be between -1 and 1. It is done, but check it out for yourself
    assertEquals(0d, annotation2.getNormalizedScore(), 0.01); // !! TODO the
scores should be between -1 and 1. It is done, but check it out for yourself
    assertEquals(-1.0, annotation3.getNormalizedScore(), 0.01); // !! !! TODO the
scores should be between -1 and 1. It is done, but check it out for yourself

}
finally {
    instance.close();
}
}
```

## 9. Running the Application

Use the following command to package and run the application if maven is in your PATH. Otherwise, use your reference framework. If IntelliJ IDEA, use the Lifecycle maven plugin to run it.

mvn clean test

### 10. Conclusion

This practice illustrates a simplified lipid annotation simulation using Drools and J2EE, demonstrating how rule engines can be integrated into scientific workflows. The project supports extensibility via new rules and metabolites and prepares the foundation for more complex reasoning and API integration.